

Est. 1989



A Mobility Shop NW Ltd.

Stairlifts - Bathrooms - Building Alterations - Scooters



Straight Stairlifts

A straightforward solution for straight staircases. The 130 stairlift is the ultimate solution to gaining your independence back around your home.

The 130 Stairlift

A straight stairlift is the ultimate staircase solution, giving you the full use of the home you love. Blending seamlessly into your home, with a smooth start and stop mechanism, you can glide effortlessly up and down the stairs.

Continually developed, improved and refined over many years, the 130 is the direct descendant of our very first straight stairlift, the Superglide. Contact us for FREE no obligation advice and speak to one of our friendly experts.

Hinged Rail Option

All our stairlifts are available with a unique hinged rail option, designed to eliminate any obstruction or trip hazard at the foot of the stairs.

New and Reconditioned Models

The 130 chairs are available to buy as either new or reconditioned models. All pre-used stairlifts are reconditioned and rigorously tested by our team of engineers, before being fitted in your home.

The Service

The reliability of our stairlifts comes from the quality of the engineering and design. For your added peace of mind, all our new and reconditioned products are backed up by a full warranty and service package. Just like any piece of machinery, we recommend that your stairlift is serviced annually by a trained engineer. We will arrange with you to come and visit to give the stairlift a thorough and complete service, putting it through its paces.



Self Charging Battery

The 130 always has the optimum amount of power so you get a smooth ride with a gentle start and stop. It even works if there is a power cut.

Swivel Seat

Using the levers enables you to swivel the seat easily, making it simple to get off the chair. Alternatively, you can choose to have a powered swivel seat, which turns the chair to face the landing automatically.

Keyswitch

Occasionally you might need to immobilise the stairlift. If children are around for example, you can use this key to make sure it can't be used.



Directional Controls

Right-handed or left-handed, we can put the controls on either arm. We can fit push-button or joystick controls, whichever you prefer.

Seatbelt

Fasten the seatbelt for added security.

Footrest Lever

To store the stairlift away, you can fold away the arm and seat. Easily creating space for those using the stairs.

Wall Control

You can use the wall controls to call the stairlift for yourself or send it to someone else.



Your Personalised Quote:

Rail Start



Inclined



Droop



Hinged



Run Off

Left Hand

Right Hand

Rail Finish



Overrun



Inclined

Extras:



Tel: 0151 220 7080
Fax: 0151 281 2115



141-145 Thomas Lane,
Broadgreen, Liverpool,
L14 5NT



www.amobilityshop.co.uk
amobshop@gmail.com