

Est. 1989



# A Mobility Shop NW Ltd.

Stairlifts - Bathrooms - Building Alterations - Scooters

## Curved Stairlifts

Designed to move safely and easily around curved and complicated turns, the 101 stairlift is compact, comfortable and convenient for awkward stairs.

### The Chair

You can choose to have push-button or joystick controls on either arm, helping you to move up and down the stairs in comfort. The swivel seat that turns to face the landing makes getting off at the top of the stairs safe and trouble-free. When not in use, the arms, seat and footrest fold away, keeping the stairs and hallway clear for other users.

### The Rail

The curved stairlifts are tailor-made to fit your individual staircase and ensure the best fit possible. The twin rails are available in a pearl colour as standard. Our rails are fitted to the stairs, not to the wall and installation is quick and clean, meaning minimal fuss and disruption in your home.

### New and Reconditioned Models

The 101 chairs are available to buy as either new or reconditioned models. All pre-used stairlifts are reconditioned and rigorously tested by our team of engineers, before being fitted in your home.

### The Service

The reliability of our stairlifts comes from the quality of the engineering and design. For your added peace of mind, all our new and reconditioned products are backed up by a full warranty and service package. Just like any piece of machinery, we recommend that your stairlift is serviced annually by a trained engineer. We will arrange with you to come and visit to give the stairlift a thorough and complete service, putting it through its paces.



### Self Charging Battery

The 101 always has the optimum amount of power so you get a smooth ride with a gentle start and stop. It even works if there is a power cut.

### Upholstery Material

You can choose from a range of colours to suit the look of your home.



### Directional Controls

Right-handed or left-handed, we can put the controls on either arm. We can fit push-button or joystick controls, whichever you prefer.



### Swivel Seat

Using the levers enables you to swivel the seat easily, making it simple to get off the chair. Alternatively, you can choose to have a powered swivel seat, which turns the chair to face the landing automatically.



### Footrest Lever

To store the stairlift away, you can fold away the arm and seat and press the lever to raise the footrest. There's no awkward bending down, which is particularly useful at the top of the stairs.



### Safety Edges

If there is something blocking the chair on the stairs, the safety edges recognise this and the lift will stop automatically.



### Wall Control

You can use the wall controls to call the stairlift for yourself or send it to someone else.



### Keyswitch

Occasionally you might need to immobilise the stairlift. If children are around for example, you can use this key to make sure it can't be used.



### Seatbelt

Fasten the seatbelt for added security.




### Curved Rails

Fitted to the stairs and not to the wall, the twin curved rails are available in a pearl colour as standard with an additional five colour options.



## Your Personalised Quote:

### Rail Start



Inclined



Droop



Hinged



Run Off

Left Hand

Right Hand

### Rail Finish



Overrun



Inclined

Extras:



Tel: 0151 220 7080  
Fax: 0151 281 2115



141-145 Thomas Lane,  
Broadgreen, Liverpool,  
L14 5NT



www.amobilityshop.co.uk  
amobshop@gmail.com